Hepalbin®-MaxiCycler

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Please read the information below before use!

English

Intended purpose/Instructions

Filter for removal of unwanted substances such as stabilizers (octanoate/N-acetyl-tryptophanate), bile acids, bilirubin and creatinine from aqueous and albumin-containing dialysate solutions.

Consultation with dialysis specialist before use is advisable.

The Hepalbin® MaxiCycler is available in different versions:

Hepalbin®-MaxiCycler 1	REF: 80005-1	1 Filter layer
Hepalbin®-MaxiCycler 3	REF: 80005-3	3 Filter layers
Hepalbin®-MaxiCycler 5	REF: 80005-5	5 Filter layers
Hepalbin®-MaxiCycler 7	REF: 80005-7	7 Filter layers
Hepalbin®-MaxiCycler 9	REF: 80005-9	9 Filter layers

Indication

Purification of albumin-containing and non-albumin-containing dialysate.

Patient group

Patients with pathological accumulation of toxins (with varying degrees of albumin binding).

Application restrictions

The use of the Hepalbin®-MaxiCycler is restricted to medical solutions that are approved as medicinal products or dialysate solution for infusion or for extracorporeal treatment.

Disposal

The Hepalbin®-MaxiCycler should be discarded after use as clinical hazardous waste in accordance with local regulations.

Contraindications/ side effects/ risks

Are the same as for extracorporeal organ support and blood purification procedures.

For further instructions see product information¹

Warnings and precautions:

- Apply careful hygienic precautions during setup.
- Should only be used up to the specified expiry date (see label).
- Designed for single use. Resterilisation or reuse leads to material failure.
- Application must be avoided or discontinued upon the detection of damage or leakage.
- Must not be exposed to internal pressures higher than 1 bar.
- Must only be used for a maximum of 24 hours.
- If the sterile barrier is accidentally opened, the product must be used immediately and must not be stored again.
- Must only be used in device combinations for which instructions for use authorized by Albutec are available.

Caution! Application only by medically trained and qualified personnel.

Caution! Only Hepalbin®-MaxiCycler with 1 or 3 layers are approved for treatment of newborns under 5 kg.

Caution! Use only in a dialysate circuit separated from the patient's blood by a dialysis membrane. Intravenous administration of dialysate without a particle filter after the Hepalbin® MaxiCycler must be excluded.

Caution! Discard the product if the sterile barrier is damaged.

Caution! risks and hazards correspond to those of other extracorporeal procedures, e. g. catheter infections.

Mandatory reporting

Serious incidents in connection with the product must be reported to Albutec GmbH and the responsible authority immediately.

For further information see product information¹

Symbol glossary

 Medical device	 Production date	 Manufacturer
 Reference number	 Lot number	 Expiry date
 Store dry	 Protect from sunlight	 Product is gamma-sterilized
 Do not use if packaging is damaged	 Do not re-sterilize	 CE marking and coding for notified body
 For single use only for max. 24 hours	 Caution - the application of the product requires the attention of the operator. For this, please follow the instructions for use!	

Please read the information below before use!

English

Application setup

The Hepalbin®-MaxiCycler is to be positioned as shown in the application sketch (outlet on the side at the top!)

If required, the following adapters can be used:

For connection to Luer-Lock connectors, use enclosed adapters. (REF AD 825.1 & AD 825.2)

Prevention of air bubbles

If the Hepalbin® MaxiCycler is integrated into a system where increased outgassing is known to occur, the use of an air trap upstream of the Hepalbin®-MaxiCycler is recommended.

Application procedure

I Pre-rinsing the Hepalbin®-MaxiCycler

The system must be filled with 0.9 % saline solution or dialysate. The filling should be as bubble-free as possible (e. g. use of an air trap). The necessary amount of rinsing solution can be found in the table (see Table, I). The escaping rinsing solution must be discarded.

II Filling the Hepalbin®-MaxiCycler


Fill the system with the dialysate used (aqueous/albuminous). The filling should be as bubble-free as possible (e. g. use of an air trap). The recommended volume can be found in the table (see Table, II). The rinsing solution displaced from the system by the dialysate filling must be discarded.

III Starting the treatment

The maximum flow rate can be found in the table.

Tabular overview of the performance features:

Stage	Hepalbin®-MaxiCycler version		1	3	5	7	9
I - III	Flow rate	[ml/min]	≤50	≤150	≤250	≤350	≤450
I	Pre-rinse volume	[ml]	300	550	800	1050	1300
II	MaxiCycler fill volume	[ml]	120	250	375	500	630

Schematic representation of the integration of a Hepalbin®-MaxiCycler⁷ in a recirculating dialysate circuit to deplete undesirable substances such as octanoate, bile acids, bilirubin and creatinine.